



EMS
Microscopy
Academy

Endless Possibilities ...

Kirsch
notes

Providing Solutions:
Buffers for
Electron Microscopy

Buffers for Electron Microscopy

There are seven common buffers used for biological sample preparation for transmission electron microscopy (TEM): Phosphate (Sorensen's Phosphate, Millonig's Phosphate), Cacodylate, PHEM, S-Collidine, Tris, Hepes, and PIPES.

There are three primary functions for the use of buffers in biological sample preparation:

- 1 The ability to maintain a constant pH (7.2 – 7.4 for mammalian and 6.9 for plants) during fixation (NOTE: The reaction between proteins and non-buffered aldehydes results in a drastic lowering of pH and the generation of morphologic artifacts.)
- 2 To provide a vehicle for the fixative which contributes to the necessary osmolarity of the solution, which should be slightly hypertonic
- 3 To contain compatible ionic moieties to avoid extraction or precipitation of constituent ions

Phosphate Buffers

Phosphate Buffer is an excellent buffer with several specific recipes, physiologically compatible with cells, and is non-toxic. Can be adapted to function optimally at a variety of pH levels by altering the ratios of the monobasic and the dibasic compounds. Typically, for pH 7.3, a 1:2.3 mono-basic : di-basic ratio is used, with more monobasic for lower or more dibasic for higher pH ranges. The osmolarity is sometimes augmented by the addition of sucrose or sodium chloride. The addition of sucrose can support growths when stored for many months, even at 4°C. Calcium cannot be used with phosphate buffer since it is precipitated with citrates.

Advantages

1. It is non-toxic
2. Has no special disposal requirements (it can be poured down the drain).
3. Can be custom tailored to function optimally at a variety of pH ranges.
4. Minimal pH change with temperature
5. Components and pre-prepared solutions are low in cost.

Disadvantages

1. It is incompatible with calcium ions.

Sorensen's Phosphate

Mono and dibasic phosphates with 0.18M sucrose added. The osmolarity of the 0.1M, pH 7.2 solution is 425. Shelf life is approximately 3 months.

Millonig's Phosphate

The typical combination of mono and dibasic phosphates, with the addition of 0.5% sodium chloride. This buffer is hypertonic (pH 7.4 - 0.1M - 440mosmols) and is recommended for very hydrated tissues. **NOTE:** For marine organisms use 3% NaCl.

Cacodylate Buffer

Cacodylate Buffer is a buffer that maintains pH levels very well during fixation. This buffer contains arsenic and is potentially carcinogenic. These properties require it to be collected and disposed of according to state and federal guide lines.

Advantages

1. Suitable for use with calcium, does not form precipitates
2. Long shelf life, does not support growths

Disadvantages

1. Potential carcinogen and contains arsenic, controlled disposal required
2. Expensive

PHEM Buffer

PHEM buffer has been used primarily for tissues and cell cultures being processed for immunocytochemical studies. Most antigens, especially intra cellular ones, stain better using PHEM than those processed using PBS, it also has a more limited effect on biochemical reactions and enzymes. Very good preservation of cellular ultrastructure, notably microtubules, is also an advantage of PHEM use.

Advantages

1. Very good preservation of cellular ultrastructure
2. Limited effect on biochemical reactions and enzymes

S-Collidine Buffer

A very stable buffer typically used to buffer OsO_4 providing excellent fixation. A very strong odor and the presence of pyrimidine makes it very toxic and difficult to work with.

Advantages

1. Excellent for buffering osmium tetroxide in 2nd fixation
2. Very good stability and buffering capacity – max. pH 7.4

Disadvantages

1. Very strong odor, requires fume hood
2. Toxic pyrimidine component
3. Not useful with aldehyde primary fixative

Tris, Tris Maleate Buffer

Advantages

1. Functional at high pH and temperature ranges used for antigen retrieval
2. Improves antibody accessibility for immuno histo/cyto chemistry

Disadvantages

1. Large pH change with temperature, must be adjusted at used temperature

Cell Culture Buffers

Hepes (Good) Buffer

A zwitterionic buffer used in a variety of cell culture environments to augment the bicarbonate buffer.

Advantages

1. Good buffering range between 6.0 – 8.0
2. Limited effect on biochemical reactions and enzymes

Pipes (Good) Buffer

Another common zwitterionic cell culture buffer.

Advantages

1. It minimizes lipid loss in animal and plant samples when buffering glutaraldehyde

Quick Comparison Table for Buffers

| EMS Catalog # Product Description | Unit | Advantages | Disadvantages |
|--|--------|--|---|
| Buffers (crystalline) | | | |
| 21190 Phosphate, Mono basic | 500 g | <ul style="list-style-type: none"> Non-toxic No special disposal Can be custom tailored at a variety of pH ranges Minimal pH change with temperature Low cost components/pre-prepared solutions | <ul style="list-style-type: none"> Incompatible with calcium ions. |
| 21180 Phosphate, Di basic | 500 g | | |
| 12300 Cacodylic Acid Sodium Salt | 100 g | <ul style="list-style-type: none"> Can use with calcium, does not form precipitates Long shelf life, does not support growths | <ul style="list-style-type: none"> Potential carcinogen Contains arsenic Controlled disposal required Expensive |
| 16782 Hepes, Good | 100 g | <ul style="list-style-type: none"> Buffering range 6.0 – 8.0 Limited effect on biochemical reactions and enzymes | |
| 19240 Pipes, Good | 100 g | <ul style="list-style-type: none"> Minimizes lipid loss in animal and plant samples when buffering glutaraldehyde | |
| Buffers: Prepared Solutions | | | |
| 19340-72 Sodium Phosphate, 0.1 M, pH 7.2 | 1L | <ul style="list-style-type: none"> Non-toxic No special disposal Can be custom tailored at a variety of pH ranges Minimal pH change with temp Low cost components/pre-prepared solutions | <ul style="list-style-type: none"> Incompatible with calcium ions. |
| 11582-10 Millonig's, 0.2M | 1L | | |
| 11600-10 Sorensen's 0.2M, pH 7.2 | 1L | | |
| 11653 Sodium Cacodylate, 0.2 M, pH 7.2 | 1L | <ul style="list-style-type: none"> Can use with calcium, does not form precipitates Long shelf life, does not support growths | <ul style="list-style-type: none"> Potential carcinogen Contains arsenic Controlled disposal required Expensive |
| 11163 PHEM | 500 ml | <ul style="list-style-type: none"> Very good preservation of cellular ultrastructure Limited effect on biochemical reactions and enzymes | |
| 11520 S-Collidine | 100 ml | <ul style="list-style-type: none"> Excellent for buffering osmium tetroxide in 2nd fixation Very good stability/buffering capacity, max. pH 7.4 | <ul style="list-style-type: none"> Very strong odor, requires fume hood Toxic pyrimidine component Not useful with aldehyde primary fixative |
| 11500 S-Collidine, pH 7.4 Kit | 200 ml | | |
| 11730-06 Tris 0.2M, pH 8.0 | 500 ml | <ul style="list-style-type: none"> Functional at high pH and temperature ranges used for antigen retrieval Improves antibody accessibility for immuno histo/cyto chemistry | <ul style="list-style-type: none"> Large pH change with temperature, must be adjusted at used temperature |
| 11740 Tris-Maleate, 0.2M, pH 6.4-8.0 | 500 ml | | |
| 11494 Hepes, 0.2M, pH 7.0 – 8.0 | 500 ml | <ul style="list-style-type: none"> Buffering range 6.0 – 8.0 Limited effect on biochemical reactions and enzymes | |
| 11610 Pipes 0.3M | 500 ml | <ul style="list-style-type: none"> Minimizes lipid loss in animal and plant samples when buffering glutaraldehyde | |