

Evaporation-controlled Automated Embedding and Polymerization

Jan LM Leunissen -  Aurion

Stacie Kirsch - **E**LECTRON **M**ICROSCOPY **S**CIENCES

Al Coritz - **E**LECTRON **M**ICROSCOPY **S**CIENCES

Traditional embedding

Standard Aldehyde and OsO₄ fixed specimen

Replace water with solvent

Freeze sub specimen



Replace solvent with resin

Refresh resin

Polymerize resin

Time involved at least 24-48hrs

Considerations

- resins are 'sticky'
- exposure to hazardous/carcinogenic chemicals
- tedious handling # specimens (a.o. duplicates)
- cross contamination risk
- interrupts other activities
- boring, yet requires attention
- risk of specimen loss in final steps
- hazardous waste

How to improve?


Don't exchange dilute resin for more concentrated mixtures/pure resin

but instead

Concentrate the resin by evaporating solvent

Until all solvent is removed: pure resin

Set Temperature – Set Vacuum – Air/N2 Flush




Critical: Vacuum Pressure

!!At pressures below the saturated vapor pressure the solvent will boil!!


Therefore at any given temperature:

the vacuum pressure should NEVER be significantly lower than the saturated vapor pressure at that temperature

Vapor Pressure of Common Solvents



Saturated Vapor Pressure of Acetone



While we're at it...can we...

- (i) automate embedding from the first solvent/resin step up to pure resin, and polymerization,
- (ii) reduce solvent and resin use with no waste, except what is left in the mixing vials,
- (iii) minimize exposure to hazardous and irritating chemicals,
- (iv) reduce hands-on time to the bare minimum,
- (v) facilitate the processing of a sufficiently large number of samples in one instrument run and
- (vi) prevent cross-contamination/specimen loss.

EMS **poly** III

Evaporation-Controlled
Automated Embedding
and
Polymerization

EMS **poly** III

Evaporation-Controlled Automated Embedding and Polymerization

3 programs for embedding (solvent parameters)

2 programs for polymerization (resin parameters)

Fully user programmable

All programs:

5 steps with both global

and individual parameter settings

EMS **poly** III

Evaporation-Controlled Automated Embedding

Overview of Program Steps

- Step 1 (**S1**): removal of gas bubbles, moderate vacuum, short time
- Step 2 (**S2**): removal of bulk of solvent,
when done auto-proceed to next step
- Step 3 (**S3**): removal of trace amounts of solvent,
using a fixed number (10) of flushing cycles
with fresh air/N₂ gas
- Steps 4 (**S4**) and 5: as Step 3, especially for less volatile solvents

EMS **poly** III


Evaporation-Controlled Automated Embedding

Parameters

- Global:
- user programmable pressure triggers
 - for gas atmosphere cleanup by venting ($\Delta P1$)
 - for pumping ON and diluted vapor removal ($\Delta P2$)
 - next Step trigger ($IT > 30min$) is factory set
 - flushing cycles with fresh air/N2 gas, factory set
- Individual:
- temperature, vacuum pressure, time

EMS poly III

Pressure vs Time Diagram



EMS **poly** III

In detail Step 2

Step 2 (**S2**): remove bulk of solvent, then auto-proceed to next step

Set temperature (T)

Pump down to pressure P just above the saturated vapor pressure (no boiling!)

Solvent evaporates → Pressure increases

When pressure increase exceeds a preset value $\Delta P1$

(indicated by **red dotted line** in previous slide)

Air is admitted (solvent gas is diluted) until a preset value $\Delta P2$

(indicated by **blue dotted line** in previous slide)

And the pump is activated until pressure is P


This procedure is repeated

until the waiting time between pump cycles exceeds 30 min

then proceed automatically to Step 3 (**S3**)

EMS poly III

In detail Step (1 and) 2



! Only the first two evaporation cycles are shown in the graph!

When IT > 30 min
Proceed to Step 3

EMS poly III

In detail Step 3 (4 and 5)

Step 3 (**S3**): removal of trace amounts of solvent, then auto-proceed to next step

Temperature and Pressure setting (P) as close as possible to saturated vapor pressure, or slightly below, to facilitate evaporation

Setup as Step 2

Continue until the interval between pump cycles exceeds 30 min....

Start Flush cycle:

[Air admit while evacuating (flush, gas ballast) for 45 secs
with intermediate vacuum for 45 secs]


3 cycles, then allow ~60 min evaporation at chosen pressure P]

Total number of gas ballast cycles: 10

When finished auto-proceed to next Step or Finish

EMS poly III

In detail Step 3 (4 and 5)



EMS **poly** III

Polymerization

At end of embedding program 1, 2 or 3:

top up BEEM capsules with pure resin
or


transfer tissue blocks into silicone molds
with resin

or

leave as small resin volume

Polymerize using programs 4 or 5

Cider yeast processed with EMS **poly**III



1% GA/caco
1% KMnO₄/dw

Acetone dehydration

EMS **poly**III
Embedded
In epoxy resin

Cider yeast processed with EMS **poly**III (OCEM, New Zealand)



1000 nm

EMS **poly** III

- No exchanges of resin or solvent/resin mix
- Fully automated, except top-up with resin
- Minimum solvent and resin use
- Environment-conscious: no resin waste
- Minimal exposure to chemicals
- Minimum hands-on time
- Up to 52 specimens in one go
- For BEEMS and Silicon molds
- No specimen loss
- Reproducible results
- Ideal for multi-user facility